

Experiment 3 (M5)

The Horizontal Projectile Motion

1- Purpose

To examine the horizontal projectile motion in the inclined plane

2- Apparatus

Air table set, launcher, wooden block, ruler.

3- Theory

An object, launched horizontally with an initial velocity v_{0x} , moves under the influence of gravitation. This motion, called as horizontal projectile motion (Fig. 3.1), has two components: one is horizontal component and the other one is vertical component. Horizontal component obeys Newton's 1st Law of Equilibrium and vertical component obeys Newton's 2nd Law of Motion; 1st and 2nd experiments, respectively.




Fig 3.1 path of a horizontal projectile motion

The motion in the horizontal component (say x - axis) has an equation of motion as:

$$x(t) = v_{ox} \cdot t = v_{ox} \cdot t$$

$$\text{at } t = 0, \quad v_{ox} = v_o$$

Additionally, the motion in the vertical component (say y - axis) has an equation of motion as:

$$y(t) = v_{oy} \cdot t + \frac{1}{2} a \cdot t^2 = \frac{1}{2} a \cdot t^2$$

$$\text{at } t = 0, \quad v_{oy} = 0$$

If we insert t (from Eq. 3.1) into Eq. 3.2, finally we get the equation of motion in vertical component (y - axis) in terms of x (displacement in horizontal component):

$$y(x) = \frac{a}{2v_o^2} x^2$$

4- Procedure

- Operate air compressor and balance the air table.
- Locate the wooden block to the back leg of the air table and find the slope of the air table.
- Keep one of the discs (pucks) unmoved at the left (or right) lower side of the air table.
- Locate the launcher, whose angle is adjusted to zero, at the right (or left) upper side of the air table.
- Put the other disc (puck) to the launcher and practise launching. When you are ready, follow 6th step.
- Adjust spark timer (generator) to a proper frequency or period.
- Operate again air compressor and launch the disc (puck). Be careful while gathering data (path A): As soon as you launch the disc (puck) push the spark timer (generator) pedal.
- Just after you gather horizontal projectile motion data, locate the disc (puck) to the top of the air table and gather the free fall data (path B).
- Check out your data if it is proper.
- Give numbers to your data points 0,1, ...,5 starting from first data point.
- Find out x- and y- axes on your data sheet (apply parallel fitting).
- Find flight time t_f , range (travelled distance on x - axis) R . Using these two data, calculate the initial velocity of the disc (puck) v_{ox} .

- Starting from zero point show the displacements on y – axis (use dimensioning rules). Record them into table 3.1.




Fig. 3.2 Data dimensioning

- Determine the accelerations of path A and path B plotting $y - t^2$ graph, separately.
- Compare these result with teoretical one, 980 cm/s^2 .

5- Data

#	B path		A path	$t \pm \Delta t$	$t^2 \pm \Delta t^2$
	$x \pm \Delta x$	$y \pm \Delta y$	$y \pm \Delta y$		
0					
1					
2					
3					
4					
5					

Table 3.1 Horizontal projectile motion data

6. Questions

1. Analyse the motion in x - and y - axes. Which type of motion does the disc follow?